VivaCart

Versatility Meets Mobility

The versatile, all-new **VivaCart**[™] offers portability, flexibility, and seamless integration with most pumps.

Roll True.

Easy-glide 5 inch (12.7 cm) swivel
casters for a smoother ride
Directional lock casters offer straighter
linear movement
Vibration dampening technology
minimizes transport noise

Cut Cart Clutter.

Vertical mounting options for monitors & equipmentSide-mount storage for tubing

protects system lines

 Soft-close internal drawer for additional storage

Tailored Transport.

Customizable upper shelf securely accommodates most pumps
 Lower shelf designed to house
 MCH-1000(m) & extends to three locking positions
 Configurable options include added

storage, cord management, power-strip, & more



CardioQuip®

LBL-82 V0 ©2023, CardioQuip LLC

VivaCart[™] | VC-1000



Specifications

Material	SS304 Electropolished Stainless Steel
Dimensions (LxWxH) including poles	81.3 x 58.4 x198.1 cm (32 x 23 x 78 in)
Caster Diameter	12.7 cm (5 in)
Weight*	46.2 kg (102 lbs)
Approved Max. Load	Top Shelf: 18.1 kg (40 lbs)
	Bottom Shelf: 36.3 kg (80 lbs)

*Base cart weight. Weight will vary depending on configuration of accessories.

Standard VivaCart Configuration:





VC-A-H16 Extended Handle, 16 cm



VC-D-SH1 Side Mount Single Hook



Pole, 25 cm

VC-C-INF Infusion Pole w/

removable hooks



VC-E-RTS Removable Upper Shelf

Customize Your Build with these available Options:



VC-E-PGM Pump Guard Mounting Brackets (set of 2)



Strip, 20A

VC-B-VSB Side Mount Storage Side Mount Power Basket (JCAHO Compliant)



VC-D-SH3 Side Mount Triple Hook



VC-A-H04 Compact Handle, 4 cm



VC-D-O2E Side Mount Oxygen Tank Holder (fits type E cylinders)

CardioQuip